

# EQUIPMAK

“we equip the whole world”

## 3 PLY PROTECTIVE MASK

Equipmak Face Masks offer the broadest range of styles, features and levels of protection available today. The materials and donning attachments are sonically bonded, and all masks have enclosed nosepieces to assist in conforming to the contours of the face. All face masks are natural rubber latex-free.



These masks provide soft comfortable protection, with no pressure to ears. They are super breathable, fluid-resistant and the triple-layer construction provides a good level of bacterial filtration efficiency and also particle efficiency.


### Technical Specifications

Brand :	Equipmak
Manufacturer :	MAK Girişim
Country of Origin :	TURKEY
Application :	Procedure Mask
Color :	White
Material :	Polypropylene (Non-woven Fabric)
Size :	17,5 cm * 9,5 cm
Sterility :	NonSterile
Sewing Style :	Ultrasonic
User :	Adult

CE

FDA





## FEATURES

- Pleat style with ties
- Protective three-layer construction
- Appropriate for general patient care and isolation where exposure to blood and/or bodily fluids is not a risk
- Designed with extra lightweight materials to enhance breathability and provide a cooler, more comfortable mask
- Stretchable earloops provide a personalized fit and make it more convenient to use

## PACKAGE :

50 Pieces Per Box,  
40 Boxes (2000 Pieces) Per CTN

